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hofmann gmbh  
hauensteinstrasse 2a  
63776 mömbris  
germany

**Thank you and congratulations.**

We are delighted you chose the mystim Tension Lover and wish you lots of joy with this love toy of the most recent generation.

In order to enjoy your device unrestrictedly it is of utmost importance to read the instruction manual thoroughly and completely before first usage of the Tension Lover.

Love toys have a long tradition and fortunately nowadays are mostly out of the closet. However especially electro-stimulation toys are still subject to prejudice. Whatever love toy serves foremost one purpose: the enrichment of sexual live. With the Tension Lover this may be achieved in an absolutely new dimension. Dare to go new ways of never suspected joy!

The mystim Tension Lover is a medical product classified 2a and is subject to constant quality control as all our other products. Should you despite find any faults or would you like to communicate ideas or statements concerning this product, please feel free to notify us any time. We will gladly be of help.

Lots of joy and bubbling entertainment while exploring wishes your mystim-Team!

## 1. TENS - what does that mean?

Stimulation of nerves by means of electrical impulse through the skin is called TENS (trans-dermal electrical nerve stimulation). TENS devices are used for pain therapy and offer a thousand-fold proven alternative to medical drugs that depict a chemical intervention to the body and thus an extra strain.

The classical TENS therapy bases on the natural body-own system of information transportation. Positive, as well as negative sensations (as pain) reach the brain by the same nervous system. TENS uses the fact, that comfortable sensations travel on faster nerve strands than negative ones. That means, when the soft electrical impulses, that TENS-therapy sends out, reach the spinal cord, there is no more room left for the pain information; it can no longer be processed by the brain.

Moreover TENS effects that body-own happiness hormones are spilled, the so called endorphins. These two factors together effect that the sensation of pain often is noticeably eased or completely eliminated.

The mystim Tension Lover is in entity a medical product and may naturally be used for pain therapy by you. Find more information about the treatment of chronic pain beginning chapter 9.

## 2. Electrical stimulation and sexual practice

For some years electrical stimulation, short e-stimulation, is getting more and more popular. The erogenous areas are run through by more nerves than average. With the mystim Tension Lover you can stimulate each one to ecstasy. You can either work rapidly to achieve an orgasm or enjoy a soft trickling that slowly spreads over the entire body for hours. In all cases you will experience a new, intensive sexuality. While dildos and vibrators are mostly out of the closet and have conquered most bedrooms, electrical stimulation devices are still under a lot of prejudice. Some of those will be explored in the following.

### **2.1 Electric current and sex – exclusively reserved for sadomasochists?**

One of the most popular misunderstandings concerning e-stimulation is the idea that electrical nerve stimulation is uncomfortable and painful and therefore exclusively suitable for use in the SM- scene.

The mystim Tension Lover is equipped with two separately controlled channels offering the possibility to adjust impulse individually according to your taste. Stimulation may be varied from gentle/soft over intensive all the way to strong/painful.

TENS devices can be used in the SM-mode as well as for gentle satisfaction.

### **2.2 E-stimulation – is perverted, isn't it?**

Definitely not, e-stimulation is not perverted. Electrical stimulation is the most recent generation of love toys. Each new toy had to first prove itself and survive existing prejudice and ideas about what is „normal“ or „abnormal“. Just a couple of years ago the use of vibrators was seen as perverted. Nowadays they are accepted by society and a common part in the sexual practice of many couples.

E-stimulation is in comparison a yet young sex toy that can offer you lust and satisfaction in a completely new dimension.

### **2.3 Electrical stimulation – isn't it dangerous?**

Danger by electrical stimulation occurs only when conducted thoughtlessly. The Tension Lover may never be adhered in a manner that current streams through the brain or the heart as would be the case when stimulating both breast warts at the same time.

People with heart problems or pacemaker may use TENS devices only after consulting their physician. If these precautions and the warnings at chapter 9 are observed, then electrical stimulation is safe.

### **2.4 E-stimulation = orgasm upon pushing a button?**

Besides negative prejudice that electrical stimulation devices are subject to some positive exaggerations are circulating in connection with e-stimulation. For example a TENS device was an orgasm-machine with which in no time never suspected highs could be experienced by pushing the button.

Such exaggerated expectations most of the time are reduced to facts after first testing the electrical stimulation device.

When using an electrical love toy three things are most important: joy while experimenting, patience and practice. In E-stimulation „practices makes the master“ is commanded.

### **3. What does an impulse consist of?**

Your mystim Tension Lover treats you by means of electrical impulse. But how is impulse defined?

#### **3.1 Impulse width**

The width of impulse determines the duration of the single impulse. By variation of intensity and impulse width it is possible to stimulate different groups of nerves. Variations in impulse width setting are therefore depending on effect to be achieved and area used in absolutely recommendable.

Please find settings for impulse width in chapter 4.2.2.

#### **3.2 Impulse frequency**

The frequency of impulse states in Hertz (Hz) the intervals of time in which electrical impulses are released to the skin. When choosing impulse frequency the preferred mode of stimulation has to be considered.

To achieve a tickling, vibrating sensation a higher frequency setting (beginning at 80 Hz), where single impulse can not be felt, is recommended. For pulsating sensation, where single impulse as such is noticed, a lower frequency setting (less than 10 Hz) is recommended. Find more information on setting impulse frequency in chapter 4.2.3.

#### **3.3 Intensity**

Intensity controls the strength of the released impulse. Every individual senses impulse strength differently due to variation in tissue resistance, thickness of skin, a. s. f.. Intensity should be chosen thus that impulse is felt as thumping or tickling, but never uncomfortable or even painful. In general stimulation should not lead to muscle contraction.

Naturally you could raise impulse intensity until muscle contraction sets in. For this a high frequency, long width and strong intensity for visible contraction is recommended.

In no instance should application be uncomfortable.

Please find more information on intensity in chapter 4.2.4.

### **3.4 Treatment time**



You can employ your mystim Tension Lover for as long, or short a time you like. Especially at low intensity and high frequency of impulse stimulation may last for hours.

Please find settings for treatment time in chapter 4.2.5.

## 4. Manual control setting

In case you do not want to apply one of the preset programs, it is possible to set pulse width, frequency, and intensity manually. To do so please push first the MODE button repeatedly until the display shows as mode „N-normal treatment mode“.

### 4.1 Draft and arrangement of control elements

## **4.2 Setting controls and their function**

### **4.2.1 Starting the unit**

You start your mystim Tension Lover by turning the intensity dial clockwise until you hear a „click“ sound. By turning counter clockwise you can thus turn the device off again.

### **4.2.2 Setting the impulse width**

The two WIDTH buttons are for setting the impulse width. Pressing the upper button increases impulse width, with the lower one you decrease width, in 1  $\mu$ s steps. Keeping one of the buttons pressed increases or decreases impulse width in steps of 5  $\mu$ s. You can read the set impulse width in the middle of the display.

The importance of impulse width is described in chapter 3.1.

### **4.2.3 Setting impulse frequency**

The RATE buttons are for setting the impulse frequency. Pressing the upper button increases impulse frequency, with the lower one you decrease frequency, in 1 Hz steps. Keeping one of the buttons pressed increases or decreases impulse frequency in steps of 5 Hz. Set impulse frequency can be read at the bottom of the display.

The importance of impulse frequency is described in chapter 3.2.

### **4.2.4 Setting impulse intensity**

The two intensity dials are for setting the impulse intensity. By turning clockwise you increase intensity, counter clockwise you decrease intensity continuously. Always start with the lowest intensity setting and slowly turn higher until sensing a tickling. Treatment should not be sensed as uncomfortable or even painful.

The importance of impulse intensity is described in chapter 3.3.

### **4.2.5 Setting treatment duration**

when turning into operation the display first shows duration of the last treatment at the top, as this is automatically memorized. By pressing the TIME button you increase duration by 10 minutes.

## 5. Preset programs

Besides setting all functions manually you may choose one of the pre-selected programs. By Pressing the MODE button you can select one of seven preset programs. With each push you switch to the next lower program. The chosen Program starts blinking in the display, left.

### 5.1 B, or burst mode

In the first mode, burst (B), single bursts of 7-10 impulses individual impulses in a fixed pattern are sent. This is a combination of common TENS and low frequency TENS. In the burst mode the frequency of stimulation is set by the device and can be controlled separately through the frequency control.

### 5.2 N, or „normal mode“

The second setting, normal mode (N) offers direct control of the setting elements of the device: intensity, pulse width, and frequency are manually set.

### 5.3 M, or modulated mode

The third setting is the modulated mode (M). through continuous change in intensity this method tries to prevent immunisation of the nerve. If intensity is increased in a period of low intensity within the cycle, intensity can be slowly increased slowly by the control element, so that you will feel this as an increase.

### 5.4 MR, or modulated pulse frequency

When applying this mode, the pulse frequency decreases continually by 50 % and then increases until the initial value is reached.

### 5.5 MRW, or modulated pulse frequency and modulated pulse width

Using the fifth mode, not only the pulse frequency, but also the pulse width are subject to increase and decrease. Oscillation occurs in the opposite manner, that means by increase of pulse frequency pulse width is decreasing, and vice versa. Intensity, pulse width, and frequency can be adjusted separately.

### **5.6 SDR, or intensity and frequency modulation**

In this sixth programme de intensity and frequency increase and decrease within a cycle of 6 seconds. Increase and decrease occur in an opposite manner.

Intensity of stimulation decreases by 20 % while frequency increases by 45 %, before returning to the initial value. Intensity, frequency , and pulse width can be adjusted separately.

### **5.7 SDW, or intensity and pulse width modulation**

The seventh programme contains a increase and decrease in intensity and pulse width within a cycle of 6 seconds. Increase and decrease occur in opposite oscillation. Stimulation intensity decreases by 20 % while pulse width increases by 54 % until they return to the initial value. Intensity, pulse width, and frequency are separately adjustable.

## **6. General safety and precautionary measures**

Do read the instruction manual, before using your mystim Tension Lover for the first time.

- Only use your TENS device according to its purpose described in the instruction manual.
- Only employ accessories recommended by mystim.
- Persons with electronic or metallic implants (f. ex. pace maker) may not be treated by TENS without prior permission of their physician.
- Do not apply the device when operating machinery or driving a vehicle.
- Please always turn off your device before application and removal of the electrodes.
- TENS devices do not have a AP/APG-safety. They may not be operated in the environment of explosive or combustible material.
- Persons afflicted with heart disease may only apply the mystim Tension Lover under extreme caution. Please consult your physician concerning side effects before application of the device.
- TENS may not be applied during pregnancy.
- The purposeful use of TENS is nerve stimulation. If the device is applied to unsuitable body areas danger to the user may arise. Never place electrodes in the laryngeal or pharyngeal area. Muscle spasms in this area can afflict asphyxiation.
- In some rare occasions skin irritations may be caused by the electrodes. In such incidence please cease stimulation and remove the electrodes. Discontinue treatment until the cause for the skin irritation has been diagnosed.

- Do not adhere electrodes to eyes, in the mouth, or over the ears.
- The electrodes are to be used only on healthy skin. Never apply electrodes to irritated or injured skin.
- Keep away from children.
- Never apply electrodes in a manner that will trans-duct currency through the brain.
- In case you want to use your mystim Tension Lover for pain therapy, you may only do so after clearing the source of pain.
- Attention: Do not plug the electrode wires into a alternating current socket!

**Additionally in Germany it is applicable as follows:**

- The device is manufactured according to the medical products law (MPG).
- The device is subject to the medical products operator's regulation (MPBetreibV).
- Operator and user are obligated to strictly observe the operator's regulation in all aspects.

For all other states of the European Community national laws and regulations for operation of medical products are valid.

## 7. Helpful hints for skin care

In order to prevent skin irritations we suggest as follows:

- Wash the areas where the electrodes are to be placed prior to application and after removal with a mild soap and water. Rinse the soap thoroughly and dry the skin.
- Excessive hair growth may be cut away with the scissors; do not shave the area that is to be stimulated.
- Many skin problems arise due to the electrodes being pulled too tight over the skin at application. To prevent this apply the electrodes from their centre to the rim and do not pull them over the skin.
- To prevent irritation by stretching tape the excess wires in a circle to the skin so it will not pull on the electrodes.
- When removing the electrodes always do so in the direction of hair growth.
- It may be helpful to apply skin lotion on the areas that electrodes were attached to between applications.
- Never adhere electrodes to skin already irritated or injured.
- Dispose of electrodes that do not adhere well.
- Self-adhesive electrodes are to be used by the same person only.
- In case of skin irritation cease usage and consult your physician.
- Prior to application read the instruction manual for self-adhesive electrodes.
- If stimulation is uncomfortable either reduce intensity to a tolerable amount or discontinue stimulation and ask your physician.

## 8. The first time

Now you are acquainted with the most important functions of your mystim Tension Lover and are ready to make your first practical experiences in the subject of electrical stimulation. Given you are already knowledgeable in e-stimulation, you most likely have your individual settings you prefer. If you are an e-stimulation new-comer, you should take the following tips as what they are: suggestions! – never meant as strict orders, as it is allowed, what pleases, and this is individually different.

### 8.1 Placemenet of electrodes

Placement of electrodes depends on many different factors: What are your personal preferences, and for you the most erogenous areas? Do you want to use your Tension Lover by yourself or with a partner? Are you male or female? Following we will make some suggestions for a bubbling experience.

#### 8.1.1 Placements for Him and Her:

Use one pair of electrodes per person. One electrode is adhered to the thigh, the other to the breast wart. Please do observe that this placement is unsuitable for persons with heart problems, and never connect both breast warts through the electrodes!

After attaching the electrodes you can set an impulse intensity, that is comfortable to both of you by slowly turning higher. Take your time and enjoy the soft trickling. You may hug and kiss and while letting the excitement rise. If you like, after some time, you may change the setting. With a little bit of sensuality you will offer your partner completely new joys.

#### 8.1.2 Placement for Him:

You can attach an electrode pair on the left and a pair on the right side of your penis. One pad each on the left and the right of the glans, and one pad each left and right on the shaft. To most e-Stimulation users the sensation is more intensive in this placement, when the currency flows from glans direction penis base.

Another option consists of placing two electrodes to one buttock, and just the remaining two at the penis.

In every instance you should take your time. Sit or lay down comfortably, and slowly turn up the impulse intensity until the vibration is most stimulating to you. When choosing a low impulse intensity the trickling sensation changes slowly to a pulsating one. You may enjoy this feeling for hours, if you like.

Many male TENS users are surprised, that their penis retreats during stimulation, although excitement is unchanged. The reason is the electrical stimulation of the erectile body of the penis. You can contravene by using a penis ring. Also „dirty minds“ or exciting pictures or movies may be of help to have the penis rise to full height.

Sometimes this mode of erectile body stimulation causes an ejaculation without orgasm.

### **8.1.3 Placement for Her.**

Especially the pelvic floor and the labia are suitable adhesion points. When stimulating the clitoris or the vagina directly, this has to be done extremely carefully, as these areas are very sensitive to pain. For first experiments a placement of one electrode on each thigh inside and one on perineum and labia is recommendable.

Lay down by yourself or with your partner and carefully turn the intensity higher, and enjoy. Additionally you can treat yourself with a vibrator or be treated, or have intercourse during electrical stimulation.

## **9. Fighting pain with the mystim Tension Lover**

Your mystim Tension Lover is a high quality medical product, that can also be used for pain therapy. In the following you will read more about correct electrode placement for a most effective treatment of pain. In case you would like to use your device more often in this area, a corresponding instruction manual can be downloaded in the net under [www.mystim.com](http://www.mystim.com).

### **9.1 Electrode placement**

Best possible electrode placement is the key to successful treatment. Besides the most usual placement strategies described as follows, it may be effective to use a deflective electrode application. Thus it is important that you take your time in experimenting for the best electrode placement for you, as the positioning is individual. Please have a TENS firm physician consult you.

#### **9.1.1 The bordering placement**

This is the most common placement strategy. Here the painful area is bordered by four electrodes. You can decide whether the current is to flow directly through the pain area or surrounds it. In order to obtain a direct flow through the spot of pain you have to place the positive terminals of the electrode pairs diagonally to the negative ones; if trying to achieve a surrounding flow, they have to be placed parallel. This is very effective if the pain is in an extremity or deep within the tissue.

#### **9.1.2 Placement in the area of a spinal nerve**

The spinal nerves exit the spinal cord in pairs between the vertebrae. In this technique of electrode placement stimulation happens above the innervated area by placing one electrode exactly on the spot of pain and the other at the nerve's root at the spinal cord.

#### **9.1.3 Placement at motive, trigger, and acupuncture points**

Motive, trigger, and acupuncture points are places of high tissue conductivity. In order to destine the exact position of such a point please consult your physician.

### **9.1.4 Multiple placement strategies**

As your TENS device possesses two independent channels, you can apply different strategies in electrode placement at the same time. One channel can be used for direct stimulation of the pain area, whereas the second channel is applied for point therapy by bordering placement.

### **10. What to observe about the electrodes?**

Do observe the instructions on electrode placement in the instruction manual in order to keep stimulation effective during treatment time and to prevent skin irritations. You may only use electrodes with the CE-certification sign.

The wires delivered with the TENS device are plugged into the sockets at the device. To do so place the device in your palm and push the plugs at the end of the wire into one of the sockets; you may use one or two connection wires. After connecting the wires please connect the electrodes to the wire ends (positive and negative terminal). Please be careful when attaching or disconnecting the wires. Do not pull at the wires, as this may lead to breakage.

### **10.1 Use of self-adhesive electrodes**

#### **10.1.1 Application**

- Clean the skin thoroughly with soap and water before application of the electrodes
- Excessive hair should be cut, but not shaved
- Connect the wires to the electrodes
- Pull the safety cover off the electrodes and apply firmly to the skin area by pressing.

#### **10.1.2 Removal**

- Lift the electrodes at the corner and gently pull off. Do not pull at the electrode wire, as this might damage the electrodes
- Re-glue the electrodes to the safety cover and then remove the wires by carefully turning and pulling.

### 10.1.3 Care and storage

Your mystim Tension Lover will only achieve best results when electrodes and wires are cared for accordingly.

- Clean the wires with a damp cloth; slight powdering with talcum (f. ex. baby powder) will prevent tangling of the wires and improves their durability
- After each Application please store the electrodes in their re-closeable bags in a cool dry place (f. ex. the refrigerator)
- Life time of the electrodes is improved by applying a drop of cool tap water (please do not use distilled water) to the adhesive area and air drying the surface turned up. Moistening the adhesive area too much may spoil the adhesive capability.

## 11. Handling the mystim Tension Lover

Your mystim Tension Lover is already equipped with a battery, but may also be operated by storage batteries.

### 11.1 Usage of re-chargeable batteries

Prior to usage of a new device the storage battery has to be charged according to the manufacturer's prescription. Please read all Instructions and safety precautions printed on the accumulator and in this leaflet before usage. Stored for more than 60 days the re-chargeable battery may discharge. Therefore prior to usage they need to be charged again, if they were stored for a duration of time.

#### 11.1.1 Charging the accumulator

- Plug the accumulator into a 220/230V socket. Using an accumulator not delivered as accessory may cause fire, electric shock or body injury.
- Regarding charging time observe the instructions of the manufacturer.
- After elapse of charging time unplug the accumulator and remove the charged battery.
- Re-chargeable batteries should always be stored loaded.

To guarantee best battery performance the following rules should be observed:

- Overloading of up to 24 hours does not damage the re-chargeable battery. Though frequent overloading may reduce their life time.

- Do not short plug the accumulator terminals. This causes the accumulator to run hot and may cause permanent damage. Do not keep batteries in clothing or purses as terminals may contact coins, keys, or other metal objects.

### **Safety precautions:**

- Do not charge any other batteries with the accumulator than re-chargeable Nickel-Cadmium-Batteries. Other tape batteries may leak or explode.
- Do not burn re-chargeable batteries, as they may explode!

### **11.2 Testing and exchanging batteries**

Delivery already contains a 9V battery. If the „low battery“ sign shows in the LCD-display, the battery too weak to operate the device and needs to be exchanged. The device automatically turns off and can not be turned on again until a new battery has been inserted. In order to guarantee save device function from time to time the battery needs to be replaced.

When replacing the battery proceed as follows:

- Pull off the battery compartment's lid downwards.
- Remove the consumed battery from the battery compartment.
- Insert a new 9V battery according to the drawing in the bottom of the battery compartment. The destined polarity has to match.
- Push the lid back upwards.

### **11.3 Care, transportation, and storage**

- The device may be cleaned with alcohol. Attention: When handling combustible fluids do not smoke or operate with open flame (e.g. candles, etc.)
- Stains and massive soiling may be removed with a mild detergent.
- Do not immerse the device into liquids or subject it to large amounts of water.
- For safe transportation always insert device into the carrying case.
- If the device is not used for a prolonged period of time, remove the batteries from the battery compartment (acid could leak and damage the device). Store the device and all accessories cool and dry in the carrying case.

### 11.4 Safety-technical control measures, STC

According to the medical products operator's regulation, § 6 (safety technical control measures), the user is obliged to have regular safety-technical controls conducted. mystim prescribes these controls in accordance to MPBetreibV § 6 in a 24-months turn.

The safety-technical control may only be delegated to persons, who are able to conduct controls in an orderly fashion based on their education, knowledge, and experience achieved by practical work, and who also have suitable measuring and proofing tools available.

Standards for safety-technical control re to be requested from the quality management agent of mystim.

#### **Prior to each application of the mystim Tension Lover device the following controls are to be conducted:**

1. Inspect device of surface damage:
  - Deformation of casing?
  - Damaged or defect electrode sockets?
2. Inspect device for damaged control elements:
  - Imprints still legible?
3. Inspect display elements:
  - The LCD-display has to show device settings when turned on.
4. Inspect accessories for usability:
  - Electrode wires non-damaged?
  - Electrodes non-damaged (no fissures, no worn material)?

### 11.5 Malfunctions

Should malfunctions occur, when using the TENS, please check as follows:

- Is the treatment mode set correctly? If not, please correct setting.
- Is the wire connected correctly to the device? Wires are to be completely inserted into the sockets.
- Does LCD-display show the device settings? If necessary, exchange the battery.
- Is the wire damaged? In any case of damage, exchange the wire.
- Are 2 Electrodes connected to the paired wire?

In case of any other problems, return the device to the dealer.

Do not try to repair a defect device yourself! Manipulation voids the warranty! For problems with the device or accessories please contact your dealer.

## **11.6 Disposal**

Device packaging should be disposed of in the recycling bin, metal parts of the device in the metal recycling disposal. Plastic parts, electrical parts, and electronic chips are to be disposed of as electronic garbage. Batteries are problem garbage and are to be disposed of accordingly. Return consumed batteries to your local recycling plant. Do not dispose of in regular garbage. You may obtain further information on disposal of problem material at your local waste management centre.

## **12. Contents of delivery and accessories**

### **12.1 Contents of delivery**

- 1 x mystim Tension Lover unit
- 2 x lead wires
- 4 x self-adhesive electrodes
- 1 x battery, 9 V
- 1 x instruction manual
- 1 x carrying case

### **Caution:**

The above mentioned contents are standard equipment, that can vary within distributors.

### **12.2 Attachments**

- self-adhesive electrodes in different forms and sizes
  - permanent electrodes for multiple use
  - gel for use with permanent electrodes
  - electrode wire for permanent electrodes
  - electrode wire for self-adhesive electrodes
  - different dildo models for vaginal and/or anal use
  - different probes for vaginal and/or anal use
  - additional attachments can be found on our internet page at [www.mystim.com](http://www.mystim.com)
- If you would like to retrofit your device to permanent electrodes, you need the fitting electrode wires and electrode gel or gel pads.

### 13. Technical specifications

Channels:	dual, isolated between channels
Pulse Amplitude:	adjustable to maximal 80 mA for positive impulse or to about 10 mA peak for negative impulse at 500 Ohm load per channel.
Pulse Rate:	adjustable from 1 to 160 Hz in steps of 1 Hz
Pulse Width:	adjustable from 50 to 260 $\mu$ s in steps of 1 $\mu$ s.
Timer:	10, 20, 30, 40, 50 and 60 minutes
Dual flank control by software:	increase of amplitude at change of modus (anti shock)
Treatment duration counter:	shows time of device usage by patient
Keylock:	prevents manipulation of parameters set by physician through user
Wave Form:	asymmetrical bi-phasic square pulse
LCD-display:	for display of set parameters
Power Source:	9V battery
Battery Life:	approximately 60 hours at regular settings
Operating Conditions:	0 °C to +50 °C at 20 to 75% relative humidity, however with no condensation, pressure 700 hPa to 1060 hPa.
Storage Conditions:	0 °C to +50 °C at 20 to 75% relative humidity, however with no condensation, pressure 700 hPa to 1060 hPa.
Dimensions:	128 mm (height) x 70 mm (width) x 25 mm (length)
Weight:	170 g (inclusive of battery)
Safety for electrical shock, type BF.	
All electrical specifications valid at $\pm$ 20 % 500 Ohm load.	

#### 13.1. Explanation of signs on the CE label

Observe the instruction manual!

Application component of type BF. Application components are not suitable to be used for direct application at the heart.

Symbol for the marking of electric appliances and electronic devices.

## **14. Warranty**

There is a 2 year guarantee for our products, applicable from date of purchase.

Within the prescribed warranty time we will remove all defects of the device, that are based on faulty material or manufacture. Removal of defect occurs at our choice by either repair or replacement of parts or the entire device.

Exempt from warranty are: Defects caused by inappropriate use (f. ex. connection to inappropriate electrical current or voltage, connection to inappropriate sources of electricity, breakage by fall and such). Normal wear and tear that only minimally affects the operational utility of the unit is not subject to the warranty.

Warranty subsides if unauthorized intervention or repair is undertaken, or if non-original parts are used.

Warranty is effective only, if date of purchase is confirmed on the warranty card by stamp and signature of the dealership.

In case of warranty claims, please return your device completely with the warranty card to one of our authorized dealerships or send both to the nearest customer service.

No other claims are valid under the provisions of this warranty, unless our liability is based on legal statute. Claims against the dealer based on the provision of the purchase agreement do not affect this warranty.

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